




# DURAN® YOUTILITY® PROMOTION LIMITED EDITION



# SPEND \$250 ON SELECTED YOUTILITY® LABORATORY BOTTLES & RECEIVE A LIMITED EDITION GENIUS BOTTLE\*



The specially shaped recessed grips on both sides of the bottle ensure easy and safe handling. The slim shape of the YOUTILITY® bottle saves space when used in autoclaves and laboratory refrigerators. The bottles can also be stored horizontally on their side. A specially marked area fits perfectly with the self-adhesive DURAN® YOUTILITY® labels.



## FEATURES + BENEFITS

- Winner of 7 International awards
- Easy handling with specially shaped gripping zones on both sides of the bottle
- Ergonomic shape of the screw cap has comfortable grooves and ridges for easier tightening or removal
- Faster opening or closing with the shortened bottle thread (GL45 compatible)
- Space saving, slimmer design. Take up less space in autoclaves.
- Made from DURAN® Borosilicate 3.3 glass

YOUTILITY® Bottles are slimmer and easier to hold than standard DURAN® bottles



# SPEND \$250 ON THE BELOW YOUTILITY® BOTTLES TO RECEIVE A LIMITED EDITION GENIUS BOTTLE\*



## Ordering Information

Supplied with Youtility® caps and pouring rings

PRODUCT CODE	VOLUME	HEIGHT	DEPTH (OD)	PACK SIZE	PRICE	SALE
<a href="#">218812854</a>	125mL	124mm	55mm	4	Login	Login
<a href="#">218813653</a>	250mL	158mm	66mm	4	Login	Login
<a href="#">218814452</a>	500mL	193mm	78mm	4	Login	Login
<a href="#">218815457</a>	1,000mL	253mm	93mm	4	Login	Login

Add name "TELSA," "CURIE" or "NEWTON" to your eligible order.

## Accessories

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE
<a href="#">292292802</a>	Cap Screw GL45 PP Cyan	PK10	Login
<a href="#">292412808</a>	Ring Pouring GL45 PP Cyan	PK16	Login
<a href="#">292432904</a>	Tag Bottle Silicone for GL45	PK16	Login

## Terms & Conditions

\*Promotion valid for retail purchases of clear glass YOUTILITY® Bottles and does not include accessories. **One** x Genius Bottle per purchase order. While stocks last only. Prices Exclude GST. E.&O.E. Select the name of the bottle you would like on your order and we will try our best to accommodate. Bottle names: "TELSA," "CURIE" or "NEWTON"

