

Balances

ANALYTICAL | HR250AZ

Our best-selling analytical balance, the HR250AZ, utilises proven A&D sensor technology to deliver superior performance with fast, accurate measurements. It has a range of 0.1 milligram to 10 micrograms readability and is covered by an industry leading 5 year warranty.

- Fast stabilisation in 2 seconds
- Large removable breeze break
- Compact footprint
- Internal calibration
- Clear, reverse-backlit LCD
- 14 weighing units
- GLP/GMP/GCP/ISO compliant output
- 90mm Ø pan size
- Percent mode and statistical calculation function

SAVE
BIG!

5 YEAR
WARRANTY

Compact & Robust
252g x 0.1mg (0.0001g)



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
HR250AZ	Balance Analytical 252g x 0.1mg (with large breeze break)	Each	LOGIN	LOGIN

COMPACT | EJ SERIES

Accurate, high quality and stackable scales are perfect for applications where portability or easy storage are paramount. The EJ Series supports multiple weighing units, has fast stabilisation and a 3 year warranty.

- AC or battery (4 x AA)
- Backlit display
- RS232 optional
- External calibration



PRODUCT CODE	DESCRIPTION	PAN SIZE	PRICE	SALE
EJ610	610g x 0.01g	110mm Ø	LOGIN	LOGIN
EJ6100	6,100g x 0.1g	127 x 140mm	LOGIN	LOGIN

PRECISION | FZ-I SERIES

Top of the range performance at an affordable price! The FZ-i series easy one touch calibration provides high reliability and fail-safe calibration.

- Fast, 1 second stabilisation with the C-SHS weight sensor
- Complies with GLP, GMP, GCP and ISO
- Multiple weighing units
- Internal calibration
- 3 year warranty
- Supplied with a small breeze break



PRODUCT CODE	DESCRIPTION	PAN SIZE	PRICE	SALE
FZ500I	520g x 0.001g	130mm Ø	LOGIN	LOGIN
FZ5000I	5200g x 0.01g	150mm Ø	LOGIN	LOGIN

pH Meters



VALUE 20 SERIES | POCKET TESTER

The economical choice for general water testing requirements, such as pools, spas, hydroponics, drinking water, cooling towers and more. Supplied with everything you need to get started!

- Large LCD display
- Self-diagnosis
- IP67 waterproof & dustproof
- Stable reading indicator

**SAVE
29%!**

PRODUCT CODE	RANGE	RESOLUTION	CALIBRATION	PACK SIZE	PRICE	SALE
PH20	0 - 14pH	0.1pH	1 - 3 points	Each	LOGIN	LOGIN



ZENTEST® SMART | POCKET TESTER



The next-generation tester designed for seamless water test experience. Professional pH measurements for general water solutions now become easier than ever.

Powered by the ZenTest® Mobile App and connected via Bluetooth®. Save and share your measurement and calibration data right on the palm of your hand via email. Data will never get lost thanks to the cloud-based data system.

- Smart self-diagnosis
- Upgraded sensor technology
- Replaceable electrode
- Cloud based data management
- Hybrid design (can be used with or without a smart phone)

PARAMETER	pH		ORP (mV)	PACK SIZE	PRICE	SALE
PRODUCT CODE	RANGE	ACCURACY	RANGE			
PH60-Z	-2 to 16 pH	±0.01 pH; ±1 digit	±1,000mV*	Each	LOGIN	LOGIN

*optional with [ORP60-DA](#) probe



850 VALUE SERIES | HANDHELD

Low cost, reliable and portable meter with simple operation.

- Large LCD display made from acrylic and is resistant to scratches
- Quick and easy 1 to 3 points automatic calibration
- IP57 waterproof rating
- Comes with everything you need, including calibration solutions and electrodes, in a convenient carry case

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
PH850	pH Meter	Each	LOGIN	LOGIN



Centrifuges



MICRO

Ideal for quick spin downs and micro filtration.

- Dual door interlock design
- Start/Stop functions by closing/opening lid
- Supplied with 1.5mL/2mL tube rotor and PCR strip rotor

Max Speed: 7,000rpm
Max RCF: 2,680 x G
Run Time: Continuous

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.003.005	Micro Centrifuge	Each	LOGIN	LOGIN



MICRO LOW SPEED

Ideal for routine medical practice in hospitals, clinics and research institutes.

- Built-in fixed angle rotor
- Speed switch (rpm/rcf)
- LCD display
- Automatic lid-lock
- Over speed detection
- Automatic internal diagnosis
- Maintenance free DC Brushless motor
- Sound alert function

Suitable for 6 tubes:

- 1.5mL/2mL
- 5mL/7mL/10mL
- 15mL

Max Speed: 5,000rpm
Range: 300 to 5,000rpm
Max RCF: 2,350 x G

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.003.020	Micro Low Speed Centrifuge	Each	LOGIN	LOGIN



MINI

The Mini's small footprint makes it perfect for every laboratory use.

- Cast iron chassis with low gravity, provides stable and reliable running
- Door interlock
- Over speed & over temperature detection
- Automatic internal diagnosis
- Sound alert function

Suitable for 12 tubes: 1.5mL/2mL

(Also 0.2mL & 0.5mL - adapters required 400.003.065 & 400.003.066)

Max Speed: 15,000rpm
Range: 500 to 15,000rpm
Max RCF: 15,100 x G

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.003.050	Mini Centrifuge	Each	LOGIN	LOGIN



Centrifuges

MINI HIGH SPEED

- Optional five rotors available
- Aluminum alloy rotors are fully autoclavable at 121°C
- LCD display
- Dual door lid-lock
- Automatic door release when rotor stops
- Over speed & over temperature detection
- Automatic internal diagnosis
- Maintenance free DC Brushless motor
- Sound alert function
- Short run time function



Max Speed: 15,000rpm
Range: 200 to 15,000rpm
Max RCF: 21,380 x G



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.003.205	Mini High Speed Centrifuge	Each	LOGIN	LOGIN

MINI HIGH SPEED REFRIGERATED

- Optional five rotors available
- Aluminum alloy rotors are fully autoclavable at 121°C
- Pre-cooling system provides cool down prior to operation
- Powerful & efficient cooling system
- Up to 9 program storage capacity
- External condensation tank prevents water accumulation within centrifuge chamber
- Electronic magnetic door locking system
- Lid remains closed when rotor stops
- Rotor imbalance detection
- Adjustable acceleration/deceleration at 9 levels allows gentle braking at low speed for better separation
- Automatic internal diagnosis
- Maintenance free DC Brushless motor
- Over speed & over temperature detection
- Sound alert and short run time function



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.003.211	Mini High Speed Refrigerated Centrifuge	Each	LOGIN	LOGIN

**SAVE
BIG!**

Magnetic Stirrers



1.5L | SQUARE PLATE

- Stirs up to 1.5L
- Square PET plate
- Compact design with small footprint
- Resistant to corrosion & easy to clean
- Supplied with electrode holder with clamp
- Maintenance free DC brushless motor

Specifications:

- Speed range: 300 - 2,000rpm
- Dimensions: 150 x 144 x 54mm
- Only 300 grams!

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.100.002	Stirrer Magnetic 1.5L LabCo	Each	LOGIN	LOGIN



20L | CIRCULAR PLATE

- Stirs up to 20L
- Ceramic coated stainless steel plate
- Resistant to corrosion & easy to clean
- Constant speed control, even during changes in load
- High magnetic field reduces stirrer bar decoupling
- Maintenance free DC brushless motor

Specifications:

- Speed range: 100 - 1,500rpm
- Dimensions: 160 x 280 x 85mm
- Weight: 2.8kg

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.100.011	Stirrer Magnetic 20L LabCo	Each	LOGIN	LOGIN



20L | SQUARE PLATE

- Stirs up to 20L
- Large square, ceramic glass plate (184mm x 184mm)
- Resistant to corrosion & easy to clean
- High magnetic field reduces stirrer bar decoupling
- Maintenance free DC brushless motor

Specifications:

- Speed range: 100 - 1,500rpm
- Dimensions: 215 x 360 x 112mm
- Weight: 5.3kg

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.100.021	Stirrer Magnetic 20L LabCo	Each	LOGIN	LOGIN




Magnetic Hotplate Stirrers



20L | CIRCULAR PLATE

- 20L stirring capacity
- Heats up to 340°C
- Speed range 100 - 1,500rpm
- Three selectable heating modes: Fast, Standard & Stable
- HOT warning, indicating residual heat
- Maintenance free DC Brushless motor
- External PT1000 temperature sensor improves accuracy to $\pm 0.2^{\circ}\text{C}$ degrees (optional extra)
- Supplied with protective cover (400.100.208)



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.100.112	Stirrer Magnetic Hotplate 20L	Each	LOGIN	LOGIN

20L | SQUARE PLATE

- Stirs up to 20L
- Large square, ceramic glass plate 184 x 184mm
- Resistant to corrosion & easy to clean
- High magnetic field reduces stirrer bar decoupling
- Maintenance free DC brushless motor

Specifications:

- Speed range: 100 - 1,500rpm
- Dimensions: 215 x 360 x 112mm
- Weight: 5.3kg



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.100.120	Stirrer Magnetic Hotplate 20L	Each	LOGIN	LOGIN

Fixed Height Trolleys



- Easy to assemble
- 100mm casters (two with brakes)
- High strength plastic shelves 715 x 455mm
- Load capacity: 275kg



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
800.010.020	Trolley Fixed Height 2 Shelves	Each	LOGIN	LOGIN
800.010.030	Trolley Fixed Height 3 Shelves	Each	LOGIN	LOGIN

Electronic Pipette Filler

Experience confidence in control and reduce fatigue with the ergonomic and lightweight LabCo® Electronic Pipette Filler. Easily aspirate up to 25mL in a matter of 5 seconds! When fully charged, you can expect more than 8 hours of intermittent use.

- 8 selectable aspiration/dispensing speeds + gravity dispensing
- Powerful and fast operation
- Large LCD displays battery life and speed setting
- Replaceable internal filter protects against over aspiration
- High quality silicone nose cone is fully autoclavable
- Low battery alert
- Can be used while charging

SPECIFICATIONS	
Aspirating Speeds	8
Dispensing Speeds	8/Gravity
Volume Range	0.1 - 100mL
BATTERY	
Battery Type	Lithium Ion
Battery Life	8 hours
Charging Time	2 hours (approx.)

Unscrew the nose cone for filter replacement and decontamination

Separate buttons for aspiration and dispensing offers precise control

Large, easily read LCD display



**SAVE
20%!**

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
550.001.005	Electronic Pipette Filler	Each	LOGIN	LOGIN

Vortex Shaker Variable Speed

- Small, touch or continuous tube mixer
- Vacuum suction feet hold the unit to the bench
- Aluminium cast base avoids unnecessary "walking"
- Versatile unit with a range of tube accessories

Specifications:

- Speed range: 0 - 2,500rpm
- Mixing diameter: 4mm
- Dimensions: 127 x 130 x 260mm
- Weight: 3.5kg

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE	SALE
400.101.001	Shaker Vortex Variable Speed	Each	LOGIN	LOGIN



lab
co
SCIENTIFIC

Vacuum Pumps

Rocker oil-free vacuum pumps are compact in size, lightweight and oil-mist-free. Perfect for vacuum filtration, microbiological tests, air sampling and suspended solids testing. Selected models are equipped with Prote™ overflow protection and a moisture trap with filter cartridge and vacuum regulator.

ROCKER 300



The best seller in the laboratory for vacuum filtration, suspended solids testing, micro and air sampling.

Applications:




Max. Vacuum 99mBar
Max. Flow Rate 20L/min

- Vacuum Regulator ✓
- Filter Cartridge ✓
- Overflow Protection ✓
- Overheat Protection ✓



PRODUCT CODE	MODEL	PRICE	PRICE
167300-22	Rocker 300	LOGIN	LOGIN

Applications
Legend:



ROCKER 400

ROCKER
make filtration easy

As versatile as the 300, but with a higher flow rate of 34L/min.

Applications:



Max. Vacuum 99mBar
Max. Flow Rate 34L/min

- Vacuum Regulator ✓
- Filter Cartridge ✓
- Overflow Protection ✓
- Overheat Protection ✓



PRODUCT CODE	MODEL	PRICE	SALE
167400-22	Rocker 400	LOGIN	LOGIN

Dry Block Heaters



SAHARA 310 & 320

Made of one-piece aluminum coated with PTFE that allows multiple uses with aluminum blocks, armor beads or water. Accurate and stable temperature control that's suitable for all types of biological and life science laboratory applications.

- Choose from 1 or 2 block heaters
- Programmable timer & auto shut off with alarm
- Built-in memory can store up to 3 user definable heating programs
- Temperature range: RT +5°C to 150°C
- Resolution: 0.1°C
- Temperature accuracy: $\pm 0.2^{\circ}\text{C}$ (@ 37°C)
- Overheat protection



PRODUCT CODE	DESCRIPTION	PACK SIZE	PRICE
179310-22	Heater Dry Block 1 Block	Each	LOGIN
179320-22	Heater Dry Block 2 Block	Each	LOGIN

**SAVE BIG ON
CLEARANCE
ITEMS**

EXPLORE OUR RANGE HERE