

VITLAB delivers reliable liquid-handling solutions for accurate dispensing, titration and pipetting. Its products combine chemical resistance, robust design and user-friendly performance across a wide range of laboratory applications.

BOTTLE TOP DISPENSER

Genius

- Volume ranges from 0.2mL to 100mL
- Broad chemical compatibility
- Fast, precise volume adjustment
- Low-waste dispensing
- Built-in safety features
- Supplied with adapters for compatibility with various bottle sizes
- Easy to maintain and fully autoclavable



PRODUCT CODE	DESCRIPTION	VOLUME RANGE	PACK SIZE	PRICE	SALE
1625503VL	Bottle Top Dispenser Genius	0.2 - 2mL	Each	LOGIN	LOGIN
1625504VL	Bottle Top Dispenser Genius	0.5 - 5mL	Each	LOGIN	LOGIN
1625505VL	Bottle Top Dispenser Genius	1 - 10mL	Each	LOGIN	LOGIN
1625506VL	Bottle Top Dispenser Genius	2.5 - 25mL	Each	LOGIN	LOGIN
1625507VL	Bottle Top Dispenser Genius	5 - 50mL	Each	LOGIN	LOGIN
1625508VL	Bottle Top Dispenser Genius	10 - 100mL	Each	LOGIN	LOGIN

Platinum-iridium / Tantalum

- Designed for high-purity trace analysis
- Platinum-iridium model suitable for hydrogen fluoride
- Tantalum model suitable for hydrogen peroxide
- Volume range from 1 - 10mL
- Available with or without recirculation
- Broad chemical compatibility
- Fast, precise volume adjustment
- Built-in safety features
- Supplied with adapters for compatibility with various bottle sizes
- Easy to maintain and fully autoclavable



PRODUCT CODE	DESCRIPTION	RECIRCULATION	VOLUME RANGE	PACK SIZE	PRICE	SALE
1627515VL	Bottle Top Dispenser TA (platinum-iridium)	No	1 - 10mL	Each	LOGIN	LOGIN
1627525VL	Bottle Top Dispenser TA (platinum-iridium)	Yes	1 - 10mL	Each	LOGIN	LOGIN
1627535VL	Bottle Top Dispenser TA (tantalum)	No	1 - 10mL	Each	LOGIN	LOGIN
1627545VL	Bottle Top Dispenser TA (tantalum)	Yes	1 - 10mL	Each	LOGIN	LOGIN

ELECTRONIC BURETTE

- Ideal for laboratory titrations
- Precise, controlled liquid dispensing
- Battery-operated, portable design
- Simple calibration
- Easy to read digital volume display
- Innovative recirculation system for reduced waste
- Available in 25mL and 50mL options

**SAVE
16%**



PRODUCT CODE	DESCRIPTION	VOLUME	PACK SIZE	PRICE	SALE
1620506VL	Burette Electronic Titrator	25mL	Each	LOGIN	LOGIN
1620507VL	Burette Electronic Titrator	50mL	Each	LOGIN	LOGIN

ELECTRONIC PIPETTE FILLER

- Pipetting range from 0.1 - 200mL
- Lightweight, ergonomic design
- Adjustable speed with one-handed control
- Secure pipette grip and vapour protection
- 8-hour battery with charging indicator

**SAVE
15%**



PRODUCT CODE	DESCRIPTION	VOLUME RANGE	PACK SIZE	PRICE	SALE
1631500VL	Electronic Pipette Filler	0.1 - 200mL	Each	LOGIN	LOGIN

HANDHELD DISPENSER

Piccolo

- One-handed operation
- Ergonomic design for comfortable use
- No pipette tips required, reducing waste and costs
- 360° rotating discharge tube for easy positioning
- Ideal for dispensing tiny quantities of liquids
- Volume ranges from 100µL to 1,000µL



PRODUCT CODE	DESCRIPTION	VOLUME	PACK SIZE	PRICE	SALE
1610501VL	Handheld Dispenser Piccolo	100µL (Fixed)	Each	LOGIN	LOGIN
1610502VL	Handheld Dispenser Piccolo	200µL (Fixed)	Each	LOGIN	LOGIN
1610503VL	Handheld Dispenser Piccolo	250µL (Fixed)	Each	LOGIN	LOGIN
1610504VL	Handheld Dispenser Piccolo	500µL (Fixed)	Each	LOGIN	LOGIN
1610506VL	Handheld Dispenser Piccolo	1,000µL (Fixed)	Each	LOGIN	LOGIN
1611503VL	Handheld Dispenser Piccolo	100µL/250µL (Dual)	Each	LOGIN	LOGIN
1611506VL	Handheld Dispenser Piccolo	500µL/1,000µL (Dual)	Each	LOGIN	LOGIN
1611508VL	Handheld Dispenser Piccolo	1,000µL/2,000µL (Dual)	Each	LOGIN	LOGIN